


HANDYPAP AUTO CPAP

Technical Information



- Size in inch: 5.6 x 2.9 x 7.7
- Weight: 1.6 lbs
- Power supply: 100-240V AC
- Sound level (ISO 17510-1:2002): 27dBA
- Display: Symbols
- Mode: Auto-CPAP, CPAP
- Pressure range: 4-20 cm H₂O
- Altitude adjustment: Automatic
- Leak compensation: Automatic
- Heated Humidification mode: Humid Control
- Ramp time: 0-45 min, i-Ramp
- Comfort mode: CC+
- Resistance compensation: Specific Calibration
- Hands-free start (Auto On): Standard
- Mask Unplugged: Standard
- Respiratory events detection: Standard
- Tubing: 15 mm, 22 mm, Specific
- Data transfer: Mini usb-cable, SD card (Oxylink)
- Connectivity: Cellular modem (optional)
- MEDTECH connect: Fully enabled
- Oximetry: Nonin oximeter optional (OxyLink)

Environmental requirements

Temperature: 41° F to 104° F / 41° F to 95° F with humidifier.

Humidity: 10% - 95% relative humidity without condensation.

Storage: - 4° F to 140° F, ≤ 95% relative humidity without condensation.

Altitude: From sea level to 7 900 ft.

Electrical characteristics

Electrical supply: Input range 100-240V (50-60 Hz).

Maximum electrical power: 25 W (without humidifier).

93/42/EC Classification: Class IIa.

HANDYPAP AUTO CPAP

Ordering Information

HANDYPAP Accessories:
HANDYPAP

GK H₂O water chamber

MT-113902

MT-114600

MT-314600-11

Cables:

USB / mini USB cable 2m (10pts)

MT-214831-10

Power Supplies:

Export power supply

MT-415910-03

Travel Accessories:

Carrying case **HANDYPAP**

MT-214531-01

Cigarette lighter cable GK420

MT-213930-03

Battery cable ring term GK420

MT-213930-02

Tubes:

Patient circuit 1,80m Ø22 mm grey

MT-261000-04

Patient circuit 1,80m Ø15 mm grey

MT-261000-00

Hybernite ROC 2300

ROC2300

Hybernite HBT 300

HBT300

Power pal for Hybernite

POWERPAL

Filters:

HANDYPAP inlet air filter

MT-315940-01


DreamStar Analyze:

DreamStar Analyze Software

MT-215630-00

DreamStar Analyze clinical kit

MT-215630-01



HANDYPAP AUTO CPAP



For more information, please contact :



MEDTECH LIFE PVT. LTD.

(Formerly Nulife Global Medical Devices Pvt. Ltd.)

B6, Byculla Service Industries,

D.K. Marg, Byculla, Mumbai 400 027, India.

T: +91 22 23748372/73 / +9122 6657 8989

E: info@medtechlife.com W : www.medtechlife.com

CE 0459

BR0152-01 / 10-2017 IN

An ISO 9001: 2008
Certified Company

Made in France