

HANDYNEB SMART

PISTON NEBULIZER

Breathe
Easy
with
Nebulizer



Nebulizers are commonly used for the treatment of asthma, cystic fibrosis, COPD and other respiratory diseases or disorders. They use oxygen, compressed air or ultrasonic power to break up solutions and suspensions into small aerosol droplets that are inhaled from the mouthpiece of the device.

Air pollution in India caused 1.2 million deaths last year*

*Lancet Planetary Health 2018

Carry the compact Medtech's Nebulizer for uninterrupted memories, anytime, anywhere. Medtech's Nebulizer comes with a BPA free 2-part medicine bottle.

FEATURES



We go the extra mile with every Nebulizer.

So you can breathe easy, every time.

We provide

- Piston-Type Nebulizer
- Powerful, Heavy Duty Motor
- 100% Copper Wire Winding for a Long Shelf-Life
- Motor with 3 Ball Bearings for Smooth Running
- More Air Flow with a Big Fan and Chamber
- Comes with Adult and Child Mask
- 2 meter long PVC Tube
- BPA-free 2-Part Medicine Cup
- Polycarbonate Medicine Cup Top

Others may provide

- Piston-Type Nebulizer
- Light Duty Motor
- Aluminium Wire Winding
- Motor with 2 Ball Bearings
- Small Fan and Chamber
- Either Adult or Child Mask
- 1.5 meter PVC tube
- 3-Part Medicine Cup
- Polypropylene Medicine Cup Top

Please take advice of the authorized medical practitioner before use.

Note: The information provided is best to our knowledge and subject to change without notice.

Picture for illustration purposes only, actual product and accessories may vary.

Specifications are subject to change without notice.

* Source : 25th March 2014 | Geneva, World Health Organization 2012.

OUR PRODUCT RANGE

Nebulizer | BP Monitor (Digital/LED) | Steam Inhaler | Thermometer | Pulse Oxymeter
Oxygen Concentrator | Air Pressure Mattress | Auto CPAP | Needle Burner & Syringe Destroyer

MEDTECH LIFE PVT. LTD. (Formerly Nulife Global Medical Devices Pvt. Ltd.)

B6, Byculla Service Industries, D.K. Marg, Byculla, Mumbai 400 027, India.

☎ +9122 2374 8372 / 73 +9122 6657 8989 ✉ info@medtechlife.com

📞 Domestic: +91 82917 26885/ +91 82917 26884 / 83

📞 Export: +91 91369 23382 / +91 93263 56572

🌐 <https://bit.ly/2NYdqh3>

📘 <https://www.facebook.com/medtech.India>

🌐 [medtechlife](http://medtechlife.com)



ISO
13485

CE
2460