



www.medtechlife.com

# NEBULIZER HANDYNEB CLASSIC



**PATENTED** DESIGN **MADE IN INDIA**

## KEY FEATURES



Piston Type Compressor



Extensive Ventilation



Safety Fuse - Better Safety



2 Part Medicine Bottle (Integrated Baffle)



Low Noise Level - Less Than 60dBA



Light Weight - 1.35 Kg



Anti-Shock Case



3 Ball Bearings For Smoother Running

## SAHI CHUNO, MEDTECH LO.

### MEDTECH LIFE PVT. LTD.

Regd. Address: B6, Byculla Service Industries, D. K. Marg, Byculla, Mumbai - 400027, Maharashtra, India.  
LEGAL MANUFACTURER: Plot No. 05, Morai Industrial Park, S.NO. 315/5, Morai, Valsad, Gujarat-396185, India

For India	+91 82917 26884	sales@medtechlife.com
For Export	+91 91369 23382	info@medtechlife.com

