

INFRARED THERMOMETER

FOREHEAD AND EAR CANAL



Monitor your Temperature Manage illness

A temperature rise is generally caused due to an infection. Monitoring any temperature rise can help you manage illness.

Measure temperature with Medtech's Infrared Thermometer.

KEY FEATURES



Ear Mode



Power Saving Mode



2 colours backlight

BENEFITS



Quick Reading



Contact Free



Dual Mode C/F



Memory mode




Useful For -

Measures Body temperature quickly and without body contact.
Ideal for babies and old patients.

HOW TO USE

FE-03

Probe	Device name	Infrared Thermometer
Battery Cover	Model	FE-03
Forehead cap	Measurement mode	Forehead and ear measurement modes (Adjusted)
	Measurement site	Forehead and ear canal
	Power supply	d.c.3V , 1.5 V type AAA (LR 03) batteries
	Measuring range:	34.0-42.9°C (93.2°F~109.2°F)
	Measuring accuracy: (At laboratory conditions)	±0.2°C/0.4°F during 35.5°C-42.0°C (95.0°F-107.6°F); ±0.3°C/0.5°F 34.0°C -35.4°C(93.2°F-95.7°F) and 42.1°C -42.9°C(107.8°F-109.2°F)
	Clinical repeatability:	within ±0.3°C
	Resolution of display	0.1°C/0.1°F
	Operation condition	10-40°C(50°F-104°F) , Relative humidity ≤Rh95%, 70-106kPa
	Storage condition	-25-55°C(-13°F-131.0°F) , Relative humidity ≤Rh95%, 70-106kPa



Size	135.5*38*41mm
Weight	about 50g(without batteries)
High body temperature hint	≥38.0°C(100.4°F)
Grade of waterproof	IP22
Electric shock	Internally powered ME equipment
Applied part	Type BF applied part, including the whole unit
Mode of operation	Continuous operation
Battery life	≥1000 times
Software version	V1.1
Note: Not intended to be sterilized. Not for use in an OXYGEN RICH ENVIRONMENT	

*The above specifications are subject to change without prior notice.

Please take advice of the authorized medical practitioner before use.

Note: The information provided is best to our knowledge and subject to change without notice. Picture for illustration purposes only, actual product and accessories may vary. Specifications are subject to change without notice.

OUR PRODUCT RANGE

Nebulizer | BP Monitor (Digital/LED) | Steam Inhaler | | Thermometer (Standard / Infrared) | Pulse Oxymeter
Oxygen Concentrator | Air Pressure Mattress | Auto CPAP | Needle Burner & Syringe Destroyer

MEDTECH LIFE PVT. LTD. (Formerly Nulife Global Medical Devices Pvt. Ltd.)

B6, Byculla Service Industries, D.K. Marg, Byculla, Mumbai 400 027, India.

+9122 2374 8372 / 73 +9122 6657 8989 info@medtechlife.com

Domestic: +91 82917 26885/ +91 82917 26884 / 83

Export: +91 91369 23382 / +91 87676 84540 / +91 93263 56572