



HANDYPAP AUTO CPAP



Technical Information

- Size in inch: 5.6 x 2.9 x 7.7
- Weight: 1.6 lbs
- Power Supply: 100-240V AC
- Sound Level (ISO 17510-1:2002): 27dBA
- Display: Symbols
- Mode: Auto-CPAP, CPAP
- Pressure Range: 4-20 cm H₂O
- Altitude Adjustment: Automatic
- Leak Compensation: Automatic
- Heat Humidification Mode: Humid Control
- Ramp Time: 0-45 min, i-Ramp
- Comfort Mode: CC+
- Resistance compensation: Specific Calibration
- Hands-Free Start (Auto On): Standard
- Mask Unplugged: Standard
- Respiratory Events Detection: Standard
- Tubing: 15 mm, 22 mm, Specific
- Data Transfer: Mini USB Cable, SD card (Oxylink)
- Connectivity: Cellular Modern (Optional)
- MEDTECH connect: Fully enabled
- Oximetry: Nonin Oximeter Optional (OxyLink)



Environmental Requirements

- **Temperature:** 41°F to 104° F / 41° F to 95° F with humidifier.
- **Humidity:** 10% - 95% relative humidity without condensation.
- **Storage:** -4° F to 140° F, ≤ 95% relative humidity without condensation.
- **Altitude:** From sea level to 7,900 ft.



Key Features

- Small, compact & very lightweight design
- Universal Power Supply
- Optional Heated Humidifier & Heated Tube
- State-of-the-art Technologies & Algorithms to optimize user comfort and compliance
- Comfort Mode
- Resistance compensation for any type of Mask and Tube
- Intelligent Ramp
- Hands-free Start (Auto On)
- One Year device's memory
- Real Patient Flow over the last 16 hrs. of use
- Sophisticated Analysis Software presenting Summary, Report & Graph
- Optional Oximetry



Electrical Characteristics

- **Electrical Supply:** Input Range 100-240V (50-60Hz)
- **Maximum Electrical Power:** 25 W (without humidifier)
- **93/42/EC Classification:** Class IIa



Ordering Information

HANDYPAP Accessories

HANDYPAP	MT-114600
GK H ₂ O Water Chamber	MT-314600-11

Cables:

USB / Mini USB Cable 2m (10pts)	MT-214831-10
---------------------------------	--------------

Power Supplies:

Export Power Supply	MT-415910-03
---------------------	--------------

Travel Accessories:

Carrying Case HANDYPAP	MT-214531-01
Cigarette Lighter Cable Gk420	MT-213930-03
Battery Cable Ring Term Gk420	MT-213930-02

Tubes:

Patient Circuit 1,80m Ø22 mm grey	MT-261000-04
Patient Circuit 1,80m Ø15 mm grey	MT-261000-00
Hybernite ROC 2300	ROC2300
Hybernite HBT 300	HBT300
Powerpal for Hybernite	POWERPAL

Filters:

HANDYPAP Inlet Air Filter	MT-315940-01
---------------------------	--------------

DreamStar Analyze:

DreamStar Analyze Software	MT-215630-00
DreamStar Analyze Clinical Kit	MT-215630-01



Please take advice of the authorized medical practitioner before use.

Note: The information provided is best to our knowledge and subject to change without notice.
Picture for illustration purposes only, actual product and accessories may vary.
Specifications are subject to change without notice.



OUR PRODUCT RANGE

Nebulizer | BP Monitor (Digital/LED) | Steam Inhaler | Thermometer | Pulse Oximeter
Oxygen Concentrator | Air Pressure Mattress | Auto CPAP | Needle Burner & Syringe Destroyer

MEDTECH LIFE PVT. LTD. (Formerly Nulife Global Medical Devices Pvt. Ltd.)

B6, Byculla Service Industries, D.K. Marg, Byculla, Mumbai 400 027, India.

☎ +9122 2374 8372 / 73 +9122 6657 8989 ✉ info@medtechlife.com

☎ Domestic: +91 82917 26885/ +91 82917 26884 / 83

☎ Export: +91 91369 23382 / +91 87676 84540 / +91 93263 56572

 <https://bit.ly/2NYdqh3>

 <https://www.facebook.com/medtechmumbai>

 [medtechlife](https://www.youtube.com/medtechlife)



ISO **CE**
13485 **2460**