

# MEDTECH-DOTS

## NEEDLE BURNER & SYRINGE DESTROYER

Dispose syringes & needles. Avoid contamination and fraud.









Syringes across the world get recycled by ways unknown to medical practitioners. This causes infections, diseases, and deaths in the millions. From Aids to Hepatitis, millions get affected by chronic diseases because of recycled needles and syringes.

Immediate disposal and destruction of used syringes is the only way to eradicate any possibility of contamination, accident, or fraud.

**Dispose used syringes and needles with Medtech-Dots.**

### KEY FEATURES

- |  |   |   |
|--|---|---|
|  Single Hole Operation       |  High Grade, Imported Steel Cutter |  Removable Discharge Tray               |
|  Shock-Proof Fuse Protection |  Changeable Cutter                 |  Advance Needle Burning and Hub Cutting |

### BENEFITS

- |  |   |  |
|--|---|--|
|  Easy to Use |  Compact |  Durable |
|--|---|--|



#### Useful For -

Cutting syringe hubs to render it non-reusable.  
Burns needles across 19g to 31g.

#### HOW TO USE



01



02



03

Size	(L) 180mm (W) 11mm (H) 120mm
Weight	2 Kg Approx
Power	220 V AC, 60 Hz
Protection	2 ampere

Note: The information provided is best to our knowledge and subject to change without notice.  
Picture for illustration purposes only, actual product and accessories may vary.  
Specifications are subject to change without notice.



#### OUR PRODUCT RANGE

Nebulizer | BP Monitor (Digital/LED) | Steam Inhaler | Thermometer | Pulse Oximeter  
Oxygen Concentrator | Air Pressure Mattress | Auto CPAP | Needle Burner & Syringe Destroyer

#### MEDTECH LIFE PVT. LTD. (Formerly Nulife Global Medical Devices Pvt. Ltd.)

B6, Byculla Service Industries, D.K. Marg, Byculla, Mumbai 400 027, India.

+9122 2374 8372 / 73 +9122 6657 8989 info@medtechlife.com

Domestic: +91 82917 26885/ +91 82917 26884 / 83

Export: +91 91369 23382 / +91 87676 84540 / +91 93263 56572

<https://bit.ly/2NYdqh3>

<https://www.facebook.com/medtechmumbai>

[medtechlife](https://www.youtube.com/medtechlife)



ISO CE  
13485 2460